

5-carboxyRh110-Ub-Dha (*human Ub sequence, synthetic*)

UbiQ code : UbiQ-131
 Batch # : B01112015-001
 Amount : 50 ug, lyophilized powder
 Purity : ≥95% by RP-HPLC
 Mol. Weight : 8.95 kDa
 Storage : powder at -20°C; solution at -80°C. Please avoid multiple freeze/thaw cycles.

Productsheet

Background. UbiQ-131 is a new and first of its kind activity based probe for Ub E1, E2 and (HECT/RBR) E3 ligases.¹ It is based on the Ub sequence in which the C-terminal Gly76 has been replaced by a dehydroalanine (Dha) residue. The N-terminus is labeled with the green fluorescent 5-carboxyrhodamine110 dye (cRh110). It has been prepared by total chemical synthesis and is therefore well-defined in terms of dye site. UbiQ-131 is processed in a native manner by Ub E1, E2 and (HECT/RBR) E3 ligases and during this process it forms an electrophilic intermediate that can react with the active site Cys residue of the E1, E2 and (HECT/RBR) E3 enzyme, thereby creating a covalent bond (Figure 1).

Sequence

cRh110-MQIFVKLTGKTTITLEVPSDTIENVKAKIQDKEGIPPDQQLIFAGKQLEDGRTLSDYNIQKESTLHLVLRRLRG-Dha




Figure 1 - Mode of action of Ub-Dha activity based probes for E1-E2 and (HECT/RBR)-E3 enzymes.

Important: sample preparation

- dissolve the powder in as little DMSO as possible (e.g. 40 mg/mL)
- add this DMSO stock slowly to milliQ (please note the order of addition)
- buffer the aq. solution as desired
- final stocks of e.g. 0.5 mg/mL will contain 1.25 vol% DMSO.
- buffer exchange using 3 kDa spin filters or dialysis membrane allows total removal of DMSO if desired; this is however not required as in general <5 vol% DMSO is well tolerated by most enzymes.

General Experimental Conditions E1 labeling assay.

UBE1 or UBA6 (1 μ M) in 50 mM HEPES pH 8, 100 mM NaCl, 10 mM MgCl₂ and 250 μ M ATP was incubated with probe (30 μ M) at 37°C for 30 min. The reaction was quenched by the addition of reducing sample buffer and heating (90°C for 10 min).

General Experimental Conditions E2 labeling assay.

E2 enzyme (2.5 μ M) and UBE1 (0.63 μ M) in 50 mM HEPES pH 7.5, 100 mM NaCl, 5 mM MgCl₂ and 250 μ M ATP were incubated with probe (12.5 μ M) at 37°C for 30 min. The reaction was quenched by the addition of reducing sample buffer and heating (90°C for 10 min).

General Experimental Conditions HECT E3 labeling assay.

Nedd4L (2.5 μ M), UBE2D (0.5 μ M) and UBE1 (0.25 μ M) were incubated with probe (50 μ M) in 50 mM HEPES pH 7.5, 100 mM NaCl, 5 mM MgCl₂ and 250 μ M ATP at 30°C for 2h. The reaction was quenched by the addition of reducing sample buffer and heating (90°C for 10 min).

Literature. (1) (a) Mulder et al. *Nat. Chem. Biol.* **2016**, doi DOI: 10.1038/NCHEMBIO.2084. (b) MPC Mulder, F. El Oualid and H. Ovaa. Adenylation enzyme inhibitors. Application WO/2016/032332 and NL2015/050596