

脱ユビキチン化酵素 (DUB) 活性測定用試薬




UbiQ社では、一般的なアミド結合型と、イソペプチド結合型の2種類の脱ユビキチン化酵素活性測定用試薬をご用意しております

amide based DUB activity assay reagents ¹	isopeptide based fluorescence polarization (FP) DUB activity assay reagents ²
---	---

ユビキチンのC末端カルボキシル基と、レポーター分子のアミノ基がアミド結合を介して結合しており、脱ユビキチン化酵素 (DUB) によって分解されると蛍光を発します。

ユビキチン(Ub)/ユビキチン様タンパク質 (UbL) とTAMRA 標識ペプチドが、ネイティブのイソペプチド結合を介して結合しており、この状態では高分子量のために蛍光励起によって大きな偏光を示します。この試薬がDUBにより分解されると、標識色素はより小さい分子へ分解され、偏光度が小さくなります。本試薬は、DUB 活性をこのような偏光の違いによって検出できます。

Luc, Rh110, AMC read-out	class I FP	class II FP
--------------------------	------------	-------------



specifications	specifications	specifications
----------------	----------------	----------------

AMC: exc/emi - 380/460nm
Rh110Gly: exc/emi - 485/535nm
LUC: exc/emi - 400/505nm

- **native isopeptide bond**
- high wavelength read out
exc/emi - 550/590 nm
- proven in numerous industrial scale HTS campaigns

- **native isopeptide bond**
- high wavelength read out
exc/emi - 550/590 nm
- proven in numerous industrial scale HTS campaigns
- **introduces substrate context**
- **custom design possible**

reagents		reagents		reagents	
code	name	code	name	code	name
UbiQ-001	Ub-AMC	UbiQ-012	Ub-FP	UbiQ-043	Ub-Ub(1-14)-FP K6 linked
UbiQ-002	Ub-Rh110Gly	UbiQ-019	Nedd8-FP	UbiQ-044	Ub-Ub(4-17)-FP K11 linked
UbiQ-036	Ub-Luc	UbiQ-020	SUMO1-FP	UbiQ-045	Ub-Ub(20-33)-FP K27 linked
		UbiQ-021	SUMO2-FP	UbiQ-046	Ub-Ub(22-35)-FP K29 linked
		UbiQ-022	SUMO3-FP	UbiQ-047	Ub-Ub(26-39)-FP K33 linked
		UbiQ-073	ISG15-FP	UbiQ-048	Ub-Ub(41-54)-FP K48 linked
				UbiQ-049	Ub-Ub(56-69)-FP K63 linked
UbiQ-L03	DUB activity assay explorer panel (UbiQ-001, UbiQ-002, UbiQ-012, UbiQ-036)			UbiQ-029	K561 (Ub) FANCD2 (557-565)-FP
				UbiQ-030	K13 Ub-PTEN(5-21)-FP
				UbiQ-038	K119 Ub-H2AX(115-143)-FP
				UbiQ-039	K119 Ub-γH2AX (115-143)-Ser140(PO4)-FP

literature

- (a) Dang et al. *Biochemistry* **1998**, 37, 1868. (b) Mason et al. *Biochemistry* **2004**, 43, 6535. (c) Hassiepin et al. *Analytical Biochem* **2007**, 371, 201. (d) Orcutt et al. *Biochim Biophys Acta* **2012**, 1823, 2079.
 - (a) Huang and Aulabaugh *Methods in Molecular Biology* **2009**, 565, 127. (b) Tirat et al. *Analytical Biochemistry* **2005**, 343, 244. (c) Geurink and El Oualid et al. *ChemBiochem* **2012**, 13, 293. (d) Mevissen et al. *Cell* **2013**, 154, 169.
- **for a complete list of references we refer to the website or the product group overview.**

活性部位反応型DUBプローブ

脱ユビキチン化酵素 (DUB)用の活性部位反応型プローブは、多くのDUBの活性部位に存在するシステイン残基に反応するよう設計されており、ユビキチンのC末端にDUB阻害剤であるVME (vinyl-methyl ester) またはPA (propagylamide) が結合しています。(Figure 1)

UbiQ社ではVME (vinylmethyl ester) とPA (propagylamide) の2種類の反応基と、各反応基で4種類のN末端タグ (TAMRA、Cy5、biotin、HA---> Table 3) を付加したプローブ (Table 1) をご用意しております。

Figure 1




Table 1 - our range of DUB activity probes

	VME ¹		PA ²	
tag	code	name	code	name
-	UbiQ-005	Ub-VME	UbiQ-057	Ub-PA
HA	UbiQ-035	HA-Ahx-Ahx-Ub-VME	UbiQ-078	HA-Ahx-Ahx-Ub-PA
Biotin	UbiQ-054	Biotin-Ahx-Ub-VME	UbiQ-076	Biotin -Ahx -Ub-PA
TAMRA	UbiQ-050	TAMRA-Ub-VME	UbiQ-058	TAMRA-Ub-PA
Cy5	UbiQ-071	Cy5-Ub-VME	UbiQ-072	Cy5-Ub-PA
	UbiQ-L02	DUB probe explorer panel (our full range of all 10 DUB activity probes)		

Table 2 - specifications of VME and PA warhead

VME	PA
targets USP and UCH DUBs	targets OTU, USP and UCH DUBs *
irreversible DUB inhibitor	irreversible DUB inhibitor which forms a covalent linkage that can be cleaved by acid treatment

* According to a proteomics study the PA probes also target DUBs of the Machado-Josephin Domain family.²

Table 3 - specifications of N-terminal tags

TAMRA ^{1a}	Cy5 ^{1a}	Biotin ¹	HA ¹
fluorescence detection	fluorescence detection	affinity tag	affinity tag
<ul style="list-style-type: none"> fast (in-gel) -, sensitive -, and distinct read-out. no background labelling by cross reactivity as seen in immunoblots exc 550 nm, abs 590 nm 	<ul style="list-style-type: none"> fast (in-gel) -, sensitive -, and distinct read-out. no background labelling by cross reactivity as seen in immunoblots exc 625 nm, abs 670 nm 	<ul style="list-style-type: none"> allows pull-down with biotin binding tools strongest known non-covalent interaction detection by western-blotting 	<ul style="list-style-type: none"> allows pull-down with HA binding tools influenza epitope tag detection by western-blotting

literature

- 1 (a) de Jong et al. *ChemBioChem* **2012**, *13*, 2251. (b) Altun et al. *Chem Biol* **2011**, *18*, 1401.
 2 Ekkebus et al. *J Am Chem Soc* **2013**, *135*, 2867. (b) Sommer et al. *Bioorg Med Chem* **2013**, *21*, 2511.
 → for a complete list of references we refer to the product group overview document.

