

# UbiQ

targeting the ubiquitin system

**SUMO2 FP** (5-TAMRA-Lys(SUMO-2)-Gly, human sequence, C58A, semi-synthetic)

UbiQ code : UbiQ-021

Batch # : B29072013-001

Amount : 50 ug, 1 mM in DMSO (3.7 µL)

Purity : >95%

Mol. Weight : 13.5 kDa

Storage : upon arrival, powder at -20°C, solution at -80°C. Protect from light and avoid multiple freeze/thaw cycles.

## Productsheet



**Background.** UbiQ-021 is fluorescence polarization assay reagent for deSUMOylating enzymes. It is based on a 5-carboxytetramethylrhodamine (TAMRA, exc 550 nm, emi 590 nm) modified Lys-Gly sequence that is linked to (6His tagged) SUMO-2 via a native isopeptide bond with the lysine side-chain. Typical substrate concentrations range from 10–100 nM. Effective concentrations of deSUMOylating enzymes can range from 0.01-10 nM but depend on specific assay conditions and method of detection.

### sequence

MGSSHHHHHHSSGLEVLFGQPGS MADEKPKEGV K TENNDHINL KVAGQDGSVV QFKIKRHTPL  
SKLMKAY**A**ER QGLSMRQIRF RFDGQPINET DTPAQLEMED EDTIDVFQQQ TGG

\* N-terminal 6His tag, protease cleavage site: PreScission™ (LEVLFG▼GP)

\* Cys58 mutated to Ala



**Figure 1.** A: LC-MS analysis. Mobile phase A= 1% CH<sub>3</sub>CN, 0.1% formic acid in water (milliQ) and B= 1% water (milliQ) and 0.1% formic acid in CH<sub>3</sub>CN. Phenomenex Kinetex C18, (2.1×50 mm), 2.6 µM; flow rate = 0.5 mL/min, column T = 40°C. Gradient: 5-95%B over 3½ min. B: SDS-PAGE analysis (12% SDS-PAGE gel), lane 1= SUMO2, lane 2= UbiQ-021, M= SeeBlue® Plus2 (Invitrogen). CBB staining. C: Fluorescence Scan (Perkin Elmer ProXPRESS 2D, exc 560 emi 590 nm) UbiQ-021.

## sample preparation

- Example: prepare a stock in milliQ of 100 µM by diluting the DMSO stock 10× in milliQ; the DMSO conc is now 10 vol%.
- Add the aq. stock to your assay buffer to a final concentration of 100 nM (DMSO concentration is now 0.01 vol%).
- Typical concentrations for assays: 25–100 nM. See reference 5 for full experimental details.

**Literature.** (1) Tirat et al. *Analyt Biochem* **2005**, 343, 244. (2) Huang et al. *Methods in Molecular Biology* **2009**, 565, 127. (3) Levine et al. *Analyt Biochem* **1997**, 247, 83. (4) Faesen et al. *Chem Biol* **2011**, 18, 1550. (5) Geurink et al. *ChemBiochem*, **2012**, 13, 293. (6) Bingol et al. *Nature* **2014**, 510, 370.

## General protocol fluorescence polarization assay

FP assays were performed on a PerkinElmer Wallac EnVision 2100 Multilabel Reader with a 531 nm excitation filter and two 579 nm emission filters. Fluorescence intensities were measured in the S (parallel) and P (perpendicular) direction. FP values are given in mP (millipolarization) and calculated using the following formula (1):

$$Polarization (mP) = \frac{S - (G \cdot P)}{S + (G \cdot P)} \cdot 1000 \quad (1)$$

The confocal optics were adjusted with the average P and S values for TAMRA-Lys-Gly and the grating factor (G) was determined using a polarization value (L) for TAMRA of 50 mP using the following formula (2):

$$G = \frac{average\ S}{average\ P} \cdot \frac{1 - \left(\frac{L}{1000}\right)}{1 + \left(\frac{L}{1000}\right)} \quad (2)$$

The assays were performed in "non-binding surface flat bottom low flange" black 384-well plates (Corning) at room temperature in a buffer containing 20 mM Tris-HCl, pH 7.5, 5 mM DTT, 100 mM NaCl, 1 mg/mL 3-[(3-cholamidopropyl) dimethylammonio] propanesulfonic acid (CHAPS) and 0.5 mg/mL bovine gamma globulin (BGG). Each well had a volume of 20  $\mu$ L. Buffer and enzyme were predispensed and the reaction was started by the addition of substrate. Kinetic data was collected in intervals of 2.5 or 3 min. From the obtained polarization values ( $P_t$ ) the amount of processed substrate ( $S_t$ ) was calculated with to the following formula (3):

$$S_t = S_0 - S_0 \times \left[ \frac{P_t - P_{min}}{P_{max} - P_{min}} \right] \quad (3)$$

Where  $P_t$  is the polarization measured (in mP);  $P_{max}$  is the polarization of 100% unprocessed substrate (determined for every reagent at all used substrate concentrations);  $P_{min}$  is the polarization of 100% processed substrate;  $S_0$  is the amount of substrate added to the reaction.

From the obtained  $P_t$  values the values for initial velocities ( $v_i$ ) were calculated. The  $v_i$  values are used to determine the Michaelis-Menten constants ( $K_m$ ,  $V_{max}$  and  $k_{cat}$ ) by fitting the data according to formula (4) (where  $k_{cat} = V_{max}/[E]$ ). All experimental data was processed using Ms Excel and Prism 4.03 (GraphPad Software, Inc.).

$$v_i = \frac{V_{max} \times S_0}{K_m + S_0} \quad (4)$$