

# UbiQ

targeting the ubiquitin system

## Ubiquitin-aminoluciferin (*human sequence, synthetic*)

UbiQ code : UbiQ-036

Batch # : B01092013-001

Amount : 100 ug, lyophilized powder


Purity :  $\geq 95\%$  by RP-HPLC and SDS-PAGE

Mol. Weight : 8.82 kDa

Storage : upon arrival, powder at  $-20^{\circ}\text{C}$ ; solution at  $-80^{\circ}\text{C}$ . Please store in dark and avoid multiple freeze/thaw cycles.

## Productsheet

**Background.** UbiQ-036 is based on ubiquitin that is C-terminally functionalized with 6-aminoluciferin (Figure 1A). Upon cleavage by a DUB, the released aminoluciferin functions as a substrate for luciferase, allowing detection of luminescence as a read-out for DUB-activity.



**Figure 1.** A: UbiQ-036. B: LC-MS analysis UbiQ-036. Mobile phase A= 1%  $\text{CH}_3\text{CN}$ , 0.1% formic acid in water and B= 1% water and 0.1% formic acid in  $\text{CH}_3\text{CN}$ . Phenomenex Kinetex C18, (2.1 $\times$ 50 mm), 2.6  $\mu\text{m}$ ); flow rate= 0.5 mL/min, column T=  $40^{\circ}\text{C}$ . Gradient: 5-95% B over 3½ min.

### important: sample preparation

- dissolve the powder in as little DMSO as possible (e.g., 20 mg/mL, 2.27 mM)
- add this DMSO stock slowly to milliQ (please note the order of addition)
- buffer the aq. solution as desired

**Literature.** (1) (a) White et al. *J Am Chem Soc* **1966**, *88*, 2015. (b) Reddy et al. *J Am Chem Soc* **2010**, *132*, 13586. (2) Orcutt et al. *Biochim Biophys Acta* **2012**, *1823*, 2079. (3) El Oualid et al. *Angew Chem Int Ed* **2010**, *49*, 10149.