

UbiQ

targeting the ubiquitin system

K48 di-Ubiquitin VME (human sequence, synthetic)

UbiQ code : UbiQ-086
Batch # : B01112014-001
Amount : 50 ug, lyophilized powder
Purity : ≥90%*
Mol. Weight : 17.11 kDa
Storage : upon arrival, powder at –20°C, solution at –80°C. Please avoid multiple freeze/thaw cycles.

Productsheet

Background. UbiQ-086 is an activity-based probe for deubiquitinating enzymes (DUBs). It is based on K48 linked diUb where Lys48 has been replaced by a diamino butyric acid residue equipped with a VME type warhead - the Dab(VME) type of structure is a DUB reactive mimic of the native isopeptidic linked Lys(Gly) residue (Figure 1). UbiQ-086 can be used for activity profiling experiments and structural studies. Please note the native distance between the proximal and distal Ub is preserved as much as possible in UbiQ-086.




Figure 1. Mode of action diUb VME labeling of DUBs. Top= native diUb, bottom= diUb VME probe. B: SDS-PAGE analysis. 12% Bolt Bis-Tris Plus gel (Lifetechnologies) in MES buffer at 190V. Coomassie Brilliant Blue G-250 staining.

* Based on SDS-PAGE analysis there is some Ub(1-75) present in the sample but this does not interfere with labeling experiments with DUBs.

important: sample preparation

- dissolve the powder in as little DMSO as possible (e.g., 20 mg/mL)
- add this DMSO stock slowly to milliQ (please note the order of addition) and mix by vortexing
- buffer the aq. stock as desired
- For more details see (open-access) reference 1: <https://www.ncbi.nlm.nih.gov/pmc/articles/PMC4159580/>

Literature. (1) Mulder & El Oualid et al. *ChemBioChem* **2014**, *15*, 946.