

UbiQ

targeting the ubiquitin system

Ub-Dha (*human sequence, synthetic*)

UbiQ code : UbiQ-101

Batch # : B01042015-001

Amount : 50 ug, lyophilized powder

Purity : ≥95%

Mol. Weight : 8.58 kDa

Storage : upon arrival, powder at –20°C; solution at –80°C. Please avoid multiple freeze/thaw cycles.

Productsheet

Background. UbiQ-101 (Ub-Dha) is an activity-based probe for Ub E1, E2 and (HECT/RBR) E3 ligases. It is based on ubiquitin (Ub) in which the C-terminal Gly76 has been replaced by a dehydroalanine (Dha) residue. It is processed in a native manner by Ub E1, E2 and (HECT/RBR) E3 ligases and during this process it forms an electrophilic intermediate that can react with the active site Cys residue of the E1, E2 and (HECT/RBR) E3 enzyme, thereby creating a covalent bond (Figure 1C).

important: sample preparation

- dissolve the powder in as little DMSO as possible (20 - 40 mg/mL)
- add the DMSO stock to milliQ (please note the order of addition) and mix
- buffer the aq. solution as desired
- For full details please see open-access reference 1: <https://www.ncbi.nlm.nih.gov/pmc/articles/PMC5108872/>

General Experimental Conditions E1 labeling assay

UBE1 or UBA6 (1 μ M) in 50 mM HEPES pH 8, 100 mM NaCl, 10 mM MgCl₂ and 250 μ M ATP was incubated with probe (30 μ M) at 37°C for 30 min. The reaction was quenched by the addition of reducing sample buffer and heating (if required). Samples can be analyzed by SDS-PAGE using Coomassie staining.

General Experimental Conditions E2 labeling assay

E2 enzyme (2.5 μ M) and UBE1 (0.63 μ M) in 50 mM HEPES pH 7.5, 100 mM NaCl, 5 mM MgCl₂ and 250 μ M ATP were incubated with probe (12.5 μ M) at 37°C for 30 min. The reaction was quenched by the addition of reducing sample buffer and heating (90°C for 10 min). Samples can be analyzed by SDS-PAGE using Coomassie staining.

General Experimental Conditions HECT E3 labeling assay

Nedd4L (2.5 μ M), UBE2D (0.5 μ M) and UBE1 (0.25 μ M) were incubated with probe (50 μ M) in 50 mM HEPES pH 7.5, 100 mM NaCl, 5 mM MgCl₂ and 250 μ M ATP at 30°C for 2h. The reaction was quenched by the addition of reducing sample buffer and heating (90°C for 10 min). Samples can be analyzed by SDS-PAGE using Coomassie staining.

Please note optimal reaction conditions can vary between E2 and E3 enzymes.

It is advised to vary the ATP concentration from 250 μ M to 5 mM and determine which is best for your experiment.

Literature. (1) Mulder et al. *Nat Chem Biol* 2016, 12, 523.

UbiQ

targeting the ubiquitin system



A

MQIFVKTLTGKTITLEVEPSDTIENVKAKIQDKEGIPPDQQRLIFAGKQLEDGRTLSDYNIQKESTLHLVLRRLRG-Dha

B



C



D




Figure 1. A: sequence. B: SDS-PAGE analysis. 12% Bolt Bis-Tris Plus gel (Life technologies) and MES running buffer. CBB staining with Coomassie G-250. C: LC-MS analysis. Mobile phase A= 1% CH₃CN, 0.1% formic acid in milliQ and B= 1% milliQ and 0.1% formic acid in CH₃CN. XBridge BEH300 C18 5µm 4.6x100mm; column T= 40°C, flow= 0.8 mL/min. Gradient: 30–60%B over 6.5 min. D: Mode of action Ub-Dha activity-based probes.