

DUB assay explorer panel

UbiQ code: UbiQ-L03

Amount : 4 x 25 µg lyophilized powder

Purity : ≥95% (purified by RP-HPLC)

Storage : upon arrival, powder at -20°C, solution at -80°C. Please avoid multiple freeze/thaw cycles and protect from

Productsheet

Background. UbiQ-L03 is a panel of 4 deubiquitinating enzyme (DUB) activity assay reagents (25 ug each), prepared by total chemical synthesis.¹

- **Ub-AMC (UbiQ-001)** is a quenched fluorogenic DUB substrate based on the C-terminal derivatization of ubiquitin with 7-amido-4-methylcoumarin (AMC). Cleavage of the amide bond between Gly76 of ubiquitin and the AMC moiety results in an increase in fluorescence at $\lambda_{ex} / \lambda_{em} = 380 / 460$ nm.²
- **Ub-Rh110Gly (UbiQ-002)** is a quenched fluorogenic DUB substrate based on the C-terminal derivatization of ubiquitin with rhodamine110-Gly (Rh110Gly). Cleavage of the amide bond between Ub Gly76 and Rh110Gly results in an increase in fluorescence at $\lambda_{ex} / \lambda_{em} = 535 / 485$ nm.³
- **TAMRA-Lys(Ub)-Gly (UbiQ-012)** is a fluorescence polarization assay reagent based on a 5-carboxytetramethylrhodamine (TAMRA, $\lambda_{ex} / \lambda_{em} = 530 / 580$ nm) modified Lys-Gly sequence that is linked to ubiquitin via a native isopeptide bond with the lysine side-chain.⁴
- **Ub-aminoluciferin (UbiQ-036)** is a DUB substrate based on the C-terminal derivatization of ubiquitin with 6-aminoluciferin (Luc). Upon cleavage by a DUB, the released 6-aminoluciferin functions as a substrate for luciferase, allowing detection of luminescence as a read-out for DUB-activity.⁵

important - sample preparation.

- dissolve the powder in as little DMSO as possible (e.g. 25 mg/mL = 25 ug in 1 uL DMSO)
- add this DMSO stock slowly to milliQ (please note the order of addition)
- next buffer as desired (e.g. with 1M HEPES to 50 mM HEPES)
- for example, a stock of 100 uM in buffer (± 0.9 mg/mL, ± 4 vol% DMSO) is diluted 1000× to prepare a final standard assay solution of 100 nM substrate concentration


Literature. (1) El Oualid et al. *Angew Chem Int Ed* **2010**, 49, 10149. (2a) Dang et al. *Biochemistry* **1998**, 37, 1868. (2b) Mason et al. *Biochemistry* **2004**, 43, 6535. (3) Hassiepin et al. *Analyt Biochem* **2007**, 371, 201. (4a) Huang et al. *Methods in Molecular Biology* **2009**, 565, 127. (4b) Geurink et al. *ChemBiochem*, **2012**, 13, 293. (5a) White et al. *J Am Chem Soc* **1966**, 88, 2015. (5b) Reddy et al. *J Am Chem Soc* **2010**, 132, 13586. (5c) Orcutt et al. *Biochim Biophys Acta* **2012**, 1823, 2079.

Ub-AMC (human sequence, synthetic)

UbiQ code : UbiQ-001
 Batch # : B01012013-001
 Purity : ≥95% by RP-HPLC
 Amount : 25 ug, lyophilized powder
 Mol. Weight : 8.72 kDa

A

MQIFVKTLTGKTITLEVEPSDTIENVKAKIQDKEGIPPDQQRLLIFAGKQLEDGRTLSDYNIQKESTLHLVLRIRGG



B





Figure 1. UbiQ-001: A: sequence. B: LC-MS analysis. Mobile phase A= 1% CH₃CN, 0.1% formic acid in water and B= 1% water and 0.1% formic acid in CH₃CN. Phenomenex Kinetex C18, (2.1×50 mm, 2.6 μM); flow rate= 0.8 mL/min, column T= 40°C. Gradient: 5-95% over 3.5 min.

Ub-Rh110Gly (human sequence, synthetic)

UbiQ code : UbiQ-002
 Batch # : B01092013-001
 Purity : ≥95% by RP-HPLC
 Amount : 25 ug, lyophilized powder
 Mol. Weight : 8.93 kDa

A

MQIFVKTLTGKTITLEVEPSDTIENVKAKIQDKEGIPPDQQRLIFAGKQLEDGRTLSDYNIQKESTLHLVLRIRGG



B





Figure 2. UbiQ-002: A: sequence. B: LC-MS analysis. Mobile phase A= 1% CH₃CN, 0.1% formic acid in water and B= 1% water and 0.1% formic acid in CH₃CN. Phenomenex Kinetex C18, (2.1×50 mm, 2.6 μM); flow rate= 0.8 mL/min, column T= 40°C. Gradient: 5-95% over 3.5 min.

5-TAMRA-Lys(Ub)-Gly-OH (*human sequence, synthetic*)

UbiQ code : UbiQ-012
 Batch # : B24012013-001
 Purity : ≥95% by RP-HPLC
 Amount : 25 ug, lyophilized powder
 Mol. Weight: 9.16 kDa

A

MQIFVKLTITGKTTITLEVEPSDTIENVKAKIQDKEGIPPDQQRLIFAGKQLEDGRTLSDYNIQKESTLHLVLRGG



B





Figure 3. UbiQ-012: A: sequence. B: LC-MS analysis. Mobile phase A= 1% CH₃CN, 0.1% formic acid in water and B= 1% water and 0.1% formic acid in CH₃CN. Phenomenex Kinetex C18, (2.1×50 mm, 2.6 μM); flow rate= 0.8 mL/min, column T= 40°C. Gradient: 5-95% over 3.5 min.

Ubiquitin-aminoluciferin (human sequence, synthetic)

UbiQ code : UbiQ-036
 Batch # : B01092013-001
 Purity : ≥95% by RP-HPLC
 Amount : 25 ug, lyophilized powder
 Mol. Weight : 8.82 kDa

A

MQIFVKTLTGKTITLEVEPSDTIENVKAKIQDKEGIPPDQQRLIFAGKQLEDGRTLSDYNIQKESTLHLVLRIRGG



B




Figure 4. UbiQ-036: A: sequence. B: LC-MS analysis. Mobile phase A= 1% CH₃CN, 0.1% formic acid in water and B= 1% water and 0.1% formic acid in CH₃CN. Phenomenex Kinetex C18, (2.1×50 mm, 2.6 μM); flow rate= 0.8 mL/min, column T= 40°C. Gradient: 5-95% over 3.5 min.