

Ubiquitin (*human sequence, recombinant*)

UbiQ code : UbiQ-098
Batch # : B01015022-001
Amount : 10 mg, 0.25 mL of a 40 mg/mL solution (4.67 mM) in 10 mM HEPES pH 7.9, 100 mM NaCl.
Purity : ≥95%
Mol. Weight : 8.56 kDa
Storage : upon arrival, store solution at –20°C for use within 6 months or –80°C for long-term storage.
Please avoid multiple freeze/thaw cycles.

Productsheet

Background. UbiQ-098 is wild-type ubiquitin (Ub) which has been prepared by recombinant expression in *E. coli*.

sequence

MQIFVKTLTGKTITLEVEPSDTIENVKAKIQDKEGIPPDQQRLLIFAGKQLEDGRTLSDYNIQKESTLHLVLRRLGG



SDS-PAGE analysis.

(1) El Oualid et al. *Angew Chem Int Ed* **2010**, 49, 10149.